



Core | **BDS** **Bulk Delivery System**

The Core Separations bulk delivery system (BDS) is designed for users operating multiple SFC systems who require an intermediate supply of cryogenic CO₂ without the need for large-scale bulk storage typically provided by major gas suppliers. It offers a cost-effective solution for delivering high-pressure CO₂ to laboratories with continuous demand, while significantly reducing the need for manual handling of individual CO₂ cylinders. The Core Separations BDS can be configured to meet a broad range of requirements, from supplying CO₂ to analytical instruments through to supporting preparative SFC applications with significantly higher CO₂ demands.



Core | How BDS works?




The BDS

The bulk delivery system is designed to compress cryogenically stored liquid CO₂, initially supplied at 22 bar, making it ready for use in either Supercritical Fluid Extraction (SFE) or Supercritical Fluid Chromatography (SFC) systems. It operates at a maximum pressure of 80 bar, allowing CO₂ to be delivered at 74 bar to prevent common issues with phase change from liquid to gas when CO₂ is piped over long distances at lower pressures. The delivery pressure is adjustable to meet specific requirements. The system requires no active cooling and maintains the buffer tank at 80 bar, with a front-end pressure regulator supplying CO₂ to the instrument at the desired pressure.

VESSEL SIZES AVAILABLE

50L | 100L



**Power requirements**

400-460 V (32A) (main Panel)
220V - 240 V (13A) Remote Panel

**Pneumatic Air Pressure (bar/psi)**

6.9 bar / 100 psi, 1/4" compression inlet

**CO2 Inlet**

20 bar, 1/4" compression inlet

**Vent Line**

1/2" compression outlet

**Weight**

100 kg (depending on options)

**Chiller**

Optional

**PC & Monitor**

Minimum of 1.5 GHz, 16 GB RAM, 250 GB storage, Ethernet port for control panel, wired or wireless connection for Internet connectivity. Google Chrome browser. Monitor 21" minimum with 1920 x 1080 pixels resolution

High-Flow CO₂ Delivery Without Cooling

"Our bulk delivery system operates efficiently without the need for cooling on either the pump or the incoming CO₂ supply. Under these conditions, it can deliver up to 500 g/min, depending on the pump model."



Enhanced Performance with Integrated Cooling

Applying cooling to both the incoming CO₂ and the pump heads significantly enhances system performance—boosting flow rates by up to four times and enabling output of over 2 kg/min. This increased capacity supports greater delivery efficiency, making the system an ideal choice for high-demand applications where consistent performance and reliability are essential

Advantages

Reduced Manual Handling – Eliminates the frequent need to swap out CO₂ cylinders, improving safety and workflow efficiency.

Lower CO₂ Cost – Cryogenic CO₂ is generally more cost-effective than using individual cylinders.

Improved Pressure Stability – When CO₂ is delivered over long distances (e.g. from outside to inside the lab), pressure drops can cause a biphasic liquid-gas mix when using cylinders. This can disrupt compression in the SFC system until the gas is purged. The BDS mitigates this by delivering CO₂ at 74 bar, significantly reducing the likelihood of gas formation during transit.